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Typical material properties**After casting**

Yield strength 0,2 %	575 MPa
Elongation	4,25 %
Tensile strength	750 MPa
E-module	170 GPa
Density	8,4 g/cm ³
Corrosion resistance	< 200 µg/cm ²
Tarnish resistance	yes
Melting range (Solidus/Liquidus)	1370 °C / 1385 °C
CTE (25 – 500°C)	~ 14,1 x 10 ⁻⁶ K ⁻¹
Casting temperature	1485 °C
Hardness	350 HV 10/30
Max. firing temperature	ca. 980 °C
Delivery form	ø 9,5 x 11 mm
Packaging unit	2,2 lb (1000 g)
Applied norms	DIN EN ISO 22674:2016

Intended use

System Duro is a beryllium- and nickel free cobalt based bonding alloy for the production of crowns and bridges.

Indication

For the production of crowns and bridges for the ceramic veneering

Modelation

The modelation should be done with wax that fire without leaving residues under consideration of the standard rules of designs for dental technicians. The frame has to be designed in an anatomical reduced form. The wall thickness should be at a minimum of 0.4 mm to secure the flow out of the melt. Consider a sufficient connector (6-9 mm²). Avoid sharp edges and undercuts.

Sprue design

We recommend the design of the sprue with a bar. The horizontal sprue should have Ø 4-5mm, the sprue to the restoration should have Ø 3 mm. Single crowns should be directly connected with a sprue of Ø 4 mm with a length of 15-20mm. Connect the sprue on massive areas e.g. palatinal and avoid the center of the muffle. The pre-heating process of the muffle should be 900 °C.

Melting and casting

System Duro should be melted in a ceramic crucible. Please do not use graphite crucibles and no flux! Avoid the overheating of the melt. Prevent multiple casts of melt bottoms. The chemical and mechanical properties can only be guaranteed for new material! Melting with open flame (acetylene / oxygen) and inductive melting: Once the cylinders are melted and the cast shadow falls across the molten metal, before the oxide skin begins to split, start the casting.

Devesting and cleaning

Let the muffle cool down to room temperature (ca. 20°C), do not quench with water. Put the cooled muffle into water to avoid dust generation during the devesting. Sandblast the surface with 125µm of aluminium oxide with 3-4 bar, then clean with a steam cleaner.

Soldering / Laser welding

We recommend commercially available solders. **System Duro** parts should not be soldered with gold or palladium solder. **System Duro** is also ideally suitable for laser-welding.

Preparation before ceramic veneering

The minimum thickness of the prepared coping should not be less than 0.3 mm. It's recommended to sandblast the frames with minimum 125 µm of aluminium oxide with 3-4 bar and clean with steam cleaner. Oxide firing is not mandatory but can be done as an option for 5 minutes at 980 °C with vacuum (cleaning firing). The frame needs to be sandblasted with aluminium oxide 125 µm and 3-4 bar and clean with air steam.

Handling conditions / Safety

Metal dust is harmful to health. Use when grinding and sandblasting dust extraction and respirator with filter FFP3 EN149.

Contraindications and side effects

If the instructions are observed during the production processes, incompatibilities with CoCr alloys are extremely rare. In case of a proven allergy against an ingredient of this alloy, the alloy must not be used for safety reasons. In exceptional cases, electrochemically induced, local irritations have been reported. When different alloy groups are used, galvanic effects might occur. Please inform your dentist regarding the contra-indications and side effects.

Disposal Instructions

Consult the material safety data sheets or national regulations for disposal.

Storage conditions

Temperature, humidity or light has no effect on the product properties.

Quantity

Please consider the package.

Our information and recommendation are based on the state of the art in science and technology and has to be considered correct to the best of our knowledge and experience on this day. The above version shall replace any previous versions.