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Intended use

System PEEK Blank is a milling disk out of filled Polyetheretherketone (PEEK) for the production temporary and permanent dental restorations.

Not intended for the production of implants!

General guidelines for handling

This instruction for use includes important processing steps and recommendations for **System PEEK Blank**.

Indication

Full anatomical crowns and bridges (max. 2 pontics with sufficient dimensioned connectors (min. 12 mm²), anatomical reduced crowns and bridges for the composite veneering with sufficient dimensioned connectors (min. 10 - 12 mm²), telescopic crowns and bridges, secondary structural gingivaformer, abutments, splints, denture base, bar supported structures.

Frame design

The design should be done with appropriate CAD software. Please consider an anatomically reduced framework design for the veneering with light curing composite. The wall thickness should not be less than 0.5 mm. Choose a sufficient connector dimension for the anterior and posterior region (12 mm²).

Milling

Please follow the instructions and parameters of the respective manufacturer of CAM Software the CNC milling machine. **System PEEK Blank** is a thermoplast material and should be processed with water cooling. For milling without water cooling a dust extraction system has to be installed. Use only milling tools which are approved for the processing of polymers and thermoplastics.

Cutting out the frameworks from the Blank

Remove the milled frameworks with suitable cutting tools and smoothing the supports.

Preparation of the surface before veneering with composite / resin

It is recommended prior to the veneering to sandblast the surface with aluminium oxide (125µm, 3 bar) and careful cleaning by steam cleaner. It is mandatory to use an approved additional chemical conditioning for the frame.

Fixation

For insertion of permanent restorations it is approved to use composite cements. Consider sufficient mechanical retention by sandblasting (125µm, 3 bar) and a chemical conditioning (Bonding agent) of the surface to permit the adhesion of composite cement.

Handling conditions / Safety

Dust may cause irritation by inhalation and in contact with skin. During grinding and sandblasting it is recommended to consider an adequate extraction system and goggles.

Contraindications and side effects

If the instructions are observed during the production processes, incompatibilities with PEEK are extremely rare. In case of a proven allergy against an ingredient it should not be used for safety reasons. Please inform your dentist regarding the contra-indications and side effects.

Disposal Instructions

Consult the material safety data sheets or national regulations for disposal.

Storage conditions

Store dry and protected against light.

Quantity

Please consider the package.

Our information and recommendation are based on the state of the art in science and technology and has to be considered correct to the best of our knowledge and experience on this day. Above version shall replace any previous versions.