System MG

DE - Gebrauchsanweisung System MG	2
EN - Instruction for use System MG	3
FR - Mode d'emploi System MG	3
CZ - Návod k použití System MG	4
EL - Οδηγίες χρήσης System MG	5
ES - Instrucciones de uso de System MG	6
HR - Upute za uporabu proizvoda System MG	7
HU - Használati útmutató a System MG termékhez	8
IT - Istruzioni per l'uso di System MG	9
PL - Instrukcja stosowania stopu System MG	10
PT - Instruções de utilização do System MG	11
RO - Instrucțiune de utilizare System MG	12
SE - Bruksanvisning för System MG	13
DA - Brugsanvisning til System MG	14
ET - Kasutusjuhend System MG	15
JA - ご使用方法 System MG	16
LT – Naudojimo instrukcija System MG	17
LV - Lietošanas pamācība System MG	18
NL - Gebruiksaanwijzing System MG	19
SK – Návod na použitie zliatiny System MG	20
SL – Navodila za uporabo System MG	21
TR - System MG Kullanım Talimatı	22

DE - Gebrauchsanweisung System MG

PRODUKTNAME

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System MG

BEZEICHNUNG

DESCRIPTION / DESCRIPTIF

Nichtedelmetall-Dental-Modellgusslegierung auf Kobalt-Basis, Typ 5
NPM - Cobalt-based dental casting alloy, type 5 / Alliage dentaire à couler à base de Cobalt, non précieux, de type 5.

ABMESSUNG

MEASUREMENTS / DIMENSIONS

Ø 9,5 mm x 11 mm

Ø 8 mm x 15 mm

INHALT

Content / Contenu

1000 g

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
EN - Instruction for use System MG

Intended use

System MG is a medical device for the manufacturing of cast partial dentures.

Only for professional user (Dental Technician, Dentist).

The intended patient group provides for persons with partially or non-dentate jaws.

Indication

For the production of model cast works for removable restorations.

Contraindication

- All indications not listed under Indication.
- In case of known allergic reactions to any of the ingredients.

Design

The design has to be made on an investment model under consideration of the standard rules of designs for partial dentures. The modelation should be done with wax that fire without leaving residues.

Sprue design

The sprues should have a diameter of 3.5 mm with a length of 20-30 mm. Massive areas should additionally have a melt reservoir (ø 5mm). The cast hopper should be 4mm above the highest point of the partial denture wax design.

Melting and casting

System MG should be melted in a ceramic crucible. Please do not use graphite crucibles and no flux! Avoid the overheating of the melt. Prevent multiple casts of melt bottoms. The chemical and mechanical properties can only be guaranteed for new material. Melting with open flame (acetylene / oxygen) and inductive melting: Once the cylinders are melted and a so-called shadow falls across the molten metal, before the oxide skin begins to split, start the casting. **Maximum temperature for casting: 1480 °C**. The preheating temperature of the muffle is approx. 950°C.

Devesting

Let the muffle cool down to room temperature (ca. 20°C), do not quench with water. Put the cooled muffle into water to avoid dust generation during the devesting. Sandblast the surface with 250 µm of aluminium oxide with 3-4 bar. Clean the partial denture base with a steam cleaner.

Elaboration

The framework can be elaborated with standard carbide cutters or aluminium oxide stones and rubbers, look for smooth transitions.

Soldering / Laser welding

System MG can be soldered with all suitable solder. **System MG** parts should not be soldered with gold or palladium solder. **System MG** is also ideally suitable for laser-welding.

Handling conditions / Safety

Metal dust is harmful to health. Use when grinding and sandblasting dust extraction and respirator with filter FFP3 – DIN EN 149.

Residual risks and side effects

If the instructions are observed during the production processes, incompatibilities with non-precious dental alloys are extremely rare. In case of a proven allergy against an ingredient of this alloy, the alloy must not be used for safety reasons. In exceptional cases, electrochemically induced, local irritations have been reported. When different alloy groups are used, galvanic effects might occur. Please inform your dentist regarding the residual risks and side effects. Any serious incident that involves the product must be reported to the manufacturer and the competent authority in the accorded country.

Disinfection of the dental prosthesis before insertion

Workpieces from the dental laboratory must be subjected to immersion or spray disinfection before insertion into the patient's oral cavity and then rinsed under running water.

Single-use

The chemical and mechanical properties can only be guaranteed for new material.

Disposal Instructions

Please dispose of metal residues and dust in an environmentally friendly manner. Do not allow waste to enter groundwater, water or sewage systems. Contact waste exchanges for recycling. Outer packaging can be disposed of in paper waste.

Storage conditions

Temperature, humidity or light has no effect on the product properties.

Our information and recommendation are based on the state of the art in science and technology and has to be considered correct to the best of our knowledge and experience on this day. The above version shall replace any previous versions.

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Adentatec GmbH