# System PEEK-Blank

DE - Gebrauchsanweisung System PEEK-Blank.....	2
EN - Instruction for use System PEEK-Blank.....	3
FR - Mode d'emploi System PEEK-Blank .....	3
CZ - Návod k použití System PEEK-Blank.....	4
EL - Οδηγίες χρήσης System PEEK-Blank .....	5
ES - Instrucciones de uso de System PEEK-Blank .....	6
HR - Upute za uporabu proizvoda System PEEK-Blank .....	7
HU - Használati útmutató a System PEEK-Blank termékhez .....	8
IT - Istruzioni per l'uso di System PEEK-Blank .....	9
PL - Instrukcja stosowania stopu System PEEK-Blank .....	10
PT - Instruções de utilização do System PEEK-Blank.....	11
RO - Instrucțiune de utilizare System PEEK-Blank.....	12
SE - Bruksanvisning för System PEEK-Blank .....	13
DA - Brugsanvisning til System PEEK-Blank.....	14
ET - Kasutusjuhend System PEEK-Blank .....	15
JA - ご使用方法 System PEEK-Blank .....	16
LT – Naudojimo instrukcija System PEEK-Blank .....	17
LV - Lietošanas pamācība System PEEK-Blank .....	18
NL - Gebruiksaanwijzing System PEEK-Blank .....	19
SK – Návod na použitie polyméru System PEEK-Blank .....	20
SL – Navodila za uporabo System PEEK-Blank .....	21
TR - System PEEK-Blank Kullanım Talimatı.....	22

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# EN - Instruction for use System PEEK-Blank

## Intended use

**System PEEK-Blank** is a medical device for the production of temporary and permanent dental restorations.

## Only for professional user (Dental Technician, Dentist).

The intended patient group provides for persons with partially or non-dentate jaws.

## Indication

Full anatomical crowns and bridges (max. 2 pontics with sufficient dimensioned connectors (min. 10 - 12 mm<sup>2</sup>)  
Anatomical reduced crowns and bridges for the composite veneering with sufficient dimensioned connectors (min. 10 - 12 mm<sup>2</sup>)  
teleskopik crowns and bridges, secondary structures, gingivaformer, abutments, splints, denture base, bar supported structures.

## Contraindication

- All indications not listed under Indication.
- In case of known allergic reactions to any of the ingredients.

## Frame design

The design should be done with appropriate CAD software. Please consider an anatomically reduced framework design for the veneering with light curing composite. The wall thickness should not be less than 0.5 mm. Choose a sufficient connector dimension for the anterior and posterior region (10 - 12 mm<sup>2</sup>).

## Milling

Please follow the instructions and parameters of the respective manufacturer of CAM Software and the CNC milling machine. **System PEEK-Blank** is a thermoplast material and should be processed with water cooling. For milling without water cooling a dust extraction system have to be installed. Use only milling tools which are approved for the processing of polymers and thermoplastics.

## Cutting out the frameworks from the Blank

Remove the milled frameworks with suitable cutting tools and smoothing the supports.

## Preparation of the surface before veneering with composite / resin

It is recommended prior the veneering to sandblast the surface with aluminium oxide (min. 50 µm, 3-4 bar) and clean by steam cleaner. It is mandatory to use an additional chemical conditioning e.g. visio.link, bredent, Senden for the frame.

## Fixation

For insertion of permanent restorations, it is approved to use composite cements (z.B. Panavia SA Cement, Kuraray Medical Inc.). Consider a sufficient mechanical retention by sandblasting (min. 50 µm, 3-4 bar) and careful cleaning by steam cleaner. Use a chemical conditioning (e.g. visio.link, bredent, Senden) on the surface.

## Handling conditions / Safety

**Dust may cause irritation by inhalation and in contact with skin. When finishing and sandblasting, please use suction and respiratory protection mask with filter FFP3 - DIN EN 149.**

## Residual risks and side effects

If the instructions are observed during the production processes, incompatibilities with PEEK are extremely rare. In case of a proven allergy against an ingredient it should not be used for safety reasons. Please inform your dentist regarding the residual risks and side effects. Any serious incident that involves the product must be reported to the manufacturer and the competent authority in the accorded country.

## Disinfection of the dental prosthesis before insertion

Workpieces from the dental laboratory must be subjected to immersion or spray disinfection before insertion into the patient's oral cavity and then rinsed under running water.

## Single-use

Used milling disc should be not processed for further use as a medical device.

## Disposal Instructions

Please dispose residues and dust in an environmentally friendly manner. Do not allow waste to enter groundwater, water or sewage systems. Contact waste exchanges for recycling. Outer packaging can be disposed of in paper waste.

## Storage conditions

Store dry and protected against light.

*Our information and recommendation are based on the state of the art in science and technology and has to be considered correct to the best of our knowledge and experience on this day. The above version shall replace any previous versions.*

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