System Soft-Blank

DE - Gebrauchsanweisung System Soft-Blank	2
EN - Instruction for use System Soft-Blank	3
FR - Mode d'emploi System Soft-Blank.....	3
CZ - Návod k použití System Soft-Blank	4
EL - Οδηγίες χρήσης System Soft-Blank.....	5
ES - Instrucciones de uso de System Soft-Blank	6
HR - Upute za uporabu proizvoda System Soft-Blank.....	7
HU - Használati útmutató a System Soft-Blank termékhez	8
IT - Istruzioni per l'uso di System Soft-Blank.....	9
PL - Instrukcja stosowania stopu System Soft-Blank	10
PT - Instruções de utilização do System Soft-Blank	11
RO - Instrucțiune de utilizare System Soft-Blank	12
SE - Bruksanvisning för System Soft-Blank.....	13
DA – Brugsanvisning for System Soft-Blank.....	14
ET – Kasutusjuhend System Soft-Blank.....	15
LT – Naudojimo instrukcija System Soft-Blank	16
LV – Lietošanas pamācība System Soft-Blank.....	17
NL – Gebruiksaanwijzing System Soft-Blank.....	18
SK – Návod na použitie zliatiny System Soft-Blank.....	19
SL – Navodila za uporabo System Soft-Blank.....	20
TR -System Soft-Blank Kullanım Talimatı	21

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EN - Instruction for use System Soft-Blank

Intended use

System Soft-Blank is a medical device for the milling of crowns and bridges.

Only for professional user (Dental Technician, Dentist).

The intended patient group provides for persons with partially or non-dentate jaws.

Indication

- Crowns and wide-span bridges up to 16 units (max. 4 pontics in the anterior region and max. 3 pontics in the posterior region),
- Implant retained superstructures
- Implant bars.

Contraindication

- All indications not listed under Indication.
- In case of known allergic reactions to any of the ingredients.

Milling

System Soft-Blank is construed for CNC milling machines. Please follow the instructions and parameters of the respective manufacturer of CAM Software and the CNC milling machine.

CAD

The design should be done with appropriate CAD software. Please consider an anatomically reduced framework design for the veneering with ceramic. The wall thickness should not be less than 0.3 mm. Choose a sufficient connector dimension (6 - 9 mm²). Sharp edges and undercuts should be avoided.

Cutting out the frameworks from the Blank

Remove the milled frameworks with suitable cutting tools and smoothing the supports.

Soldering / Laser welding

System Soft-Blank can be soldered with all suitable solder. **System Soft-Blank** parts should not be soldered with gold or palladium solder. **System Soft-Blank** is also ideally suitable for laser-welding.

Preparation before ceramic veneering

The frameworks can be elaborated with standard carbide cutters, look for smooth transitions and avoid overlapping material. Please use the same cutter for one alloy to avoid contamination. The minimum thickness of the prepared coping should not be less than 0.3 mm. It's recommended to sandblast the frames with minimum 110 µm of Aluminium oxide with 3-4 bar and clean with steam cleaner. Oxide firing is not mandatory but can be done as an option for 5 minutes at 980°C (cleaning firing). The frame needs to be sandblasted again to remove the present oxide layer thoroughly. In the end the cleaning by steam cleaner is mandatory. If you use a ceramic bonder, please consider the instruction for use of the manufacturer.

Handling conditions / Safety

Metal dust is harmful to health. Use when grinding and sandblasting dust extraction and respirator with filter FFP3 – DIN EN 149.

Residual risks and side effects

If the instructions are observed during the production processes, incompatibilities with CoCr alloys are extremely rare. In case of a proven allergy against an ingredient of this alloy, the alloy must not be used for safety reasons. In exceptional cases, electrochemically induced, local irritations have been reported. When different alloy groups are used, galvanic effects might occur. Please inform your dentist regarding the residual risks and side effects. Any serious incident that involves the product must be reported to the manufacturer and the competent authority in the accorded country.

Disinfection of the dental prosthesis before insertion

Workpieces from the dental laboratory must be subjected to immersion or spray disinfection before insertion into the patient's oral cavity and then rinsed under running water.

Single-use

Used milling disc should be not processed for further use as a medical device.

Disposal Instructions

Please dispose of metal residues and dust in an environmentally friendly manner. Do not allow waste to enter groundwater, water or sewage systems. Contact waste exchanges for recycling. Outer packaging can be disposed of in paper waste.

Storage conditions

Temperature, humidity or light has no effect on the product properties

Our information and recommendation are based on the state of the art in science and technology and has to be considered correct to the best of our knowledge and experience on this day. The above version shall replace any previous versions.

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Adentatec GmbH