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**Cobalt based dental casting alloy, type 5**CE  
0123**Typical composition [%]**

Co	61,8
Cr	29,5
Mo	5,7
Si	0,95
Fe	0,75
C	0,6
Mn	0,55
No other elements	> 0,1

**Typical material properties****After casting**

Yield strength 0,2 %	590 MPa
Elongation	2,8 %
Tensile strength	850 MPa
E-module	155 GPa
Density	8,25 g/cm <sup>3</sup>
Corrosion resistance	< 200 µg/cm <sup>2</sup>
Tarnish resistance	ja
Melting range (Solidus/Liquidus)	1310 °C / 1380 °C
Casting temperature	1480 °C
Hardness	460 HV 10/30
Delivery form	ø 9,5 x 11 mm
Packaging unit	2,2 lb (1000 g)
Applied norms	DIN EN ISO 22674:2016

**Intended use**

**System MG** is a beryllium- and nickel-free, cobalt based alloy for the manufacture of partial dentures.

**Indication**

For the production of model cast works for removable restorations.

**Modelation**

The design has to be made on an investment model (**Adenta®-Vest PA**, Adentatec GmbH) under consideration of the standard rules of designs for partial dentures. The modelation should be done with wax that fire without leaving residues.

**Sprue design**

The sprues should have a diameter of 3.5 mm with a length of 20-30 mm. Massive areas should additional have a melt reservoir (ø 5mm). The cast hopper should be 4mm above the highest point of the partial denture wax design.

**Recommended pre-heating temperature: 900 – 950°C (Adenta®-Vest PA, Adentatec GmbH).**

**Melting and casting**

**System MG** should be melted in a ceramic crucible. **Please do not use graphite crucibles and no flux!** Avoid the overheating of the melt. Prevent multiple casts of melt bottoms. The chemical and mechanical properties can only be guaranteed for new material. Melting with open flame (acetylene / oxygen) and inductive melting: Once the cylinders are melted and a so-called shadow falls across the molten metal, before the oxide skin begins to split, start the casting.

**Devesting and elaboration**

Let the muffle cool down to room temperature (ca. 20°C), do not quench with water. Put the cooled muffle into water to avoid dust generation during the devesting. Sandblast the surface with 250 µm of aluminium oxide with 3-4 bar. Clean the partial denture base with a steam cleaner. The framework can be elaborated with standard carbide cutters or aluminium oxide stones and rubbers.

**Soldering / Laser welding**

**System MG** is approved for soldering and laser welding with commercial solders and laser wire. We recommend the solder **System SOL** and laser wire **System Flow** from Adentatec GmbH.

**Handling conditions / Safety**

**Metal dust is harmful to health. Use when grinding and sandblasting dust extraction and respirator with filter FFP3 EN149.**

**Contraindications and side effects**

If the instructions are observed during the production processes, incompatibilities with CoCr alloys are extremely rare. In case of a proven allergy against an ingredient of this alloy, the alloy must not be used for safety reasons. In exceptional cases, electrochemically induced, local irritations have been reported. When different alloy groups are used, galvanic effects might occur. Please inform your dentist regarding the contra-indications and side effects.

**Disposal Instructions**

Consult the material safety data sheets or national regulations for disposal.

**Storage conditions**

Temperature, humidity or light has no effect on the product properties.

**Quantity**

Please consider the package.

*Our information and recommendation are based on the state of the art in science and technology and has to be considered correct to the best of our knowledge and experience on this day. The above version shall replace any previous versions.*