System PEEK-Blank

DE - Gebrauchsanweisung System PEEK-Blank.....	2
EN - Instruction for use System PEEK-Blank.....	3
FR - Mode d'emploi System PEEK-Blank	3
CZ - Návod k použití System PEEK-Blank.....	4
EL - Οδηγίες χρήσης System PEEK-Blank	5
ES - Instrucciones de uso de System PEEK-Blank	6
HR - Upute za uporabu proizvoda System PEEK-Blank	7
HU - Használati útmutató a System PEEK-Blank termékhez	8
IT - Istruzioni per l'uso di System PEEK-Blank	9
PL - Instrukcja stosowania stopu System PEEK-Blank.....	10
PT - Instruções de utilização do System PEEK-Blank.....	11
RO - Instrucțiune de utilizare System PEEK-Blank.....	12
SE - Bruksanvisning för System PEEK-Blank	13
DA - Brugsanvisning til System PEEK-Blank.....	14
ET - Kasutusjuhend System PEEK-Blank	15
LT – Naudojimo instrukcija System PEEK-Blank	16
LV - Lietošanas pamācība System PEEK-Blank	17
NL - Gebruiksaanwijzing System PEEK-Blank	18
SK – Návod na použitie polyméru System PEEK-Blank	19
SL – Navodila za uporabo System PEEK-Blank	20
TR - System PEEK-Blank Kullanım Talimatı.....	21

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Stand / Status 01/2025 as

EN - Instruction for use System PEEK-Blank

Intended use
System PEEK-Blank is a medical device for the production of temporary and permanent dental restorations.

Only for professional user (Dental Technician, Dentist).

The intended patient group provides for persons with partially or non-dentate jaws.

Indication
Full anatomical crowns and bridges (max. 2 pontics with sufficient dimensioned connectors (min. 10 - 12 mm²)
Crowns and bridges for the composite veneering with sufficient dimensioned connectors (min. 10 - 12 mm²)
Removable restorations on copings, Superstructures.

Contraindication
In case of known allergic reactions to any of the ingredients.

Frame design
The design should be done with appropriate CAD software. Please consider an anatomically reduced framework design for the veneering with light curing composite. The wall thickness should not be less than 0.5 mm. Choose a sufficient connector dimension for the anterior and posterior region (10 - 12 mm²).

Milling
Please follow the instructions and parameters of the respective manufacturer of CAM Software and the CNC milling machine. **System PEEK-Blank** is a thermoplast material and should be processed with water cooling. For milling without water cooling a dust extraction system have to be installed. Use only milling tools which are approved for the processing of polymers and thermoplastics.

Cutting out the frameworks from the Blank
Remove the milled frameworks with suitable cutting tools and smoothing the supports.

Preparation of the surface before veneering with composite / resin
It's recommended prior the veneering to sandblast the surface with aluminium oxide (min. 50 µm, 3-4 bar) and clean by steam cleaner. It is mandatory to use an additional chemical conditioning e.g. visio.link, bredent, Senden for the frame.

Fixation
For insertion of permanent restorations, it's approved to use composite cements (z.B. Panavia SA Cement, Kuraray Medical Inc.). Consider a sufficient mechanical retention by sandblasting (min. 50 µm, 3-4 bar) and careful cleaning by steam cleaner. Use a chemical conditioning (e.g. visio.link, bredent, Senden) on the surface.

Handling conditions / Safety
Dust may cause irritation by inhalation and in contact with skin. When finishing and sandblasting, please use suction and respiratory protection mask with filter FFP3 - DIN EN 149.

Residual risks and side effects
If the instructions are observed during the production processes, incompatibilities with PEEK are extremely rare. In case of a proven allergy against an ingredient it should not be used for safety reasons. Please inform your dentist regarding the residual risks and side effects. Any serious incident that involves the product must be reported to the manufacturer and the competent authority in the accorded country. The SSCP is available at <https://ec.europa.eu/tools/eudamed>.

Disinfection of the dental prosthesis before insertion
Workpieces from the dental laboratory must be subjected to immersion or spray disinfection before insertion into the patient's oral cavity and then rinsed under running water.

Single-use
Used milling disc should be not processed for further use as a medical device.

Disposal Instructions
Please dispose residues and dust in an environmentally friendly manner. Do not allow waste to enter groundwater, water or sewage systems. Contact waste exchanges for recycling. Outer packaging can be disposed of in paper waste.

Storage conditions
Store dry and protected against light.

Our information and recommendation are based on the state of the art in science and technology and has to be considered correct to the best of our knowledge and experience on this day. The above version shall replace any previous versions.

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Adentatec GmbH