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## Phosphate bonded investment powder for crowns and bridges

**Adenta®-Vest CB** is a fine corn investment powder for universal applications, for casting alloys as well as for pressable ceramics.

**Adenta®-Vest CB** is approved for the speed heating and for the conventional heating process.

### Processing

**Mixing ratio:** 100 g Powder: 24 - 26 ml liquid  
160 g Powder: 39 - 41 ml liquid

### Concentration of

| <b>mixing liquid:</b> | Material      | Inlay,<br>Veneers | Crowns | Bridges |
|-----------------------|---------------|-------------------|--------|---------|
|                       | CoCr-alloy    | 85 %*             | 85 %*  | 95 %*   |
|                       | Gold-alloy    | 55 %*             | 60 %*  | 65 %*   |
|                       | Press ceramic | 60 %*             | 60 %*  | -       |

\*The liquid concentration is only for orientation. The final liquid concentration must be determined by means of test casts.

**Storage- and processing temperature:** ca. 20-23 °C

**Pre-heating temperature:** 850 - 900°C

### Mixing

Rinse the mixing bowl and dry with a clean towel. Add the required quantity of liquid into the mixing bowl. Add the appropriate quantity of **Adenta®-Vest CB** powder. Mix well with a spatula (30 s), then mixing under vacuum for 1 min, hold the vacuum for further 10 s after mixing.

Higher room temperatures shorten the processing time of the investment material! Note: Increase the mixing liquid by 1 ml.

### Investing

It's recommended doing the investing with silicone ring ((Use muffle fleece or roughen surfaces!)). Once an investing with metal ring would be made it's mandatory to use sufficient layer (1-2) of fleece inlays. The fill in to the casting ring should be made with the lowest power of the vibrator. Once the ring has been filled up the vibrator have to switch off. Processing time ca. 5 min.

### Speed heating

Ca. 25 minutes after investing the casting form can be put into the hot furnace with a temperature of 850-900°C. Do not open the furnace door within the first 20 min.! Risk of injury!

| <b>Ring size</b> | <b>Time [min]</b> |
|------------------|-------------------|
| 1                | 35                |
| 3                | 45                |
| 6                | 60                |

### Conventional heating

With the conventional heating, the full setting expansion will be reached, you should therefore, reduce the liquid concentration by 5% - 10%.

Put the set up casting form into the cold oven. Please keep in mind to set a holding time at 290°C and 580°C for each 30-45 min (holding time depends on ring size). Hold the final temperature (850°C-900°C) also for 30 - 45 min.

**General heating rate approx. 3 - 5 °C/min with linear oven control.**

### Safety notes

- Do not open the furnace door during heat-up.
- Investment powder contains silicone. Do not inhale the dust! Risk of lung damage (silicosis, lung cancer). Recommendation: Use a Respirator, type FFP3-EN149. Avoid dust during opening and during fill-in to the mixing bowl.
- Remove the dust at your working place only with wet towels.
- Avoid dust during deinvesting, place the cooled-down muffle in water until it is thoroughly moistened.
- Use an extraction system during sandblasting.

**Warning**

Our information and recommendation are based on the state of the art in science and technology and has to be considered correct to the best of our knowledge and experience on this day. The above version shall replace any previous versions.