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**Investment material for model casting****Processing**

**Mixing ratio:** Silicone-duplication:  
100 g powder : 20-22 ml fluid → 400 g : 80 ml

Gel-duplication:  
100 g powder : 18-20 ml fluid → 400 g : 72 ml

**Concentration of mixing liquid:**

Partial denture with brackets 50-65%\*  
Tertiary restorations: 65-75%\*  
Secondary frames: 85-100%\*

\*The liquid concentration is only for orientation. The final liquid concentration must be determined by means of test casts.

**Storage- and processing temperature:** 19 °C

**Mixing advice** 20 s manual  
50 s vacuum mixer

**Setting time** ca. 9 min

**Pre-heating temperature** 950°C for speed cast

**Mixing**

Rinse the mixing bowl and dry with a clean towel. Add the required quantity of liquid into the mixing bowl. Add the appropriate quantity of **Adenta® -Vest PA** powder. Mix well with a spatula (for 20 s), then mix under vacuum for 1 min (about 350 rpm), hold the vacuum for further 10 s. after mixing.

Hot temperatures reduce the processing time! Once the temperatures in the lab increase over 30°C, we recommend increasing the quantity of liquid e.g. for silicone duplication from 20 to 22 ml / 100g and adding 10 s to the mixing time.

**Investing**

The fill-in into the casting ring should be made at the lowest power vibrator. Once the casting ring has been filled up, the vibrator must be switched off. Processing time approx. 5 min, setting time approx. 9 min.

**The over bedding can be done with distilled water.**

**Speed heating**

25 - 30 min after fill-in, the casting ring can be put into the hot furnace at a temperature of 900 - 1000°C. Do not open the furnace door within the first 20 minutes! Risk of injury!

**Conventional heating**

With the conventional heating, the full setting expansion will be reached, you should therefore reduce the liquid concentration by 5% - 10%.

Put the set-up casting form into the cold oven. Please keep in mind to set a holding time at 280°C and 580°C of 45 - 60 min (depending on the casting ring size and number of rings). Hold the final temperature also for 30 - 45 min (depending on the alloy).

**General heating rate approx. 3 - 7 °C/min with linear oven control.**

**Safety notes**

- Do not open the furnace door during heat-up.
- Investment powder contains silicone. Do not inhale the dust! Risk of lung damage (silicosis, lung cancer).

Recommendation: Use a Respirator, type FFP3-EN149. Avoid dust during opening and during fill-in to the mixing bowl.

- Remove the dust at your working place only with wet towels.
- Avoid dust during devesting, place the cooled-down muffle in water until it is thoroughly moistened.
- Use an extraction system during sandblasting.

**Storage conditions**

Store investment powder always cool and dry.

**Quantity**

Please consider the package.

*Our information and recommendation are based on the state of the art in science and technology and has to be considered correct to the best of our knowledge and experience on this day. The above version shall replace any previous versions.*